Table of Contents

01. Safety Instructions 1

02. Device Description 3

03. Preparing for Use 4

04. How to Use 5

05. Product Management 8

06. After-Sales Service 9

Safety Instructions TB-2386

※ Important Notice: This manual is available in two languages: English and Croatian. If any dispute or confusion arises caused by the translation, the English version shall prevail.

※ To ensure the proper operation of this device, please read this instruction manual carefully before using the device and keep it for further reference.

⚠️ WARNINGS

The following persons should not use this device:

- Persons with wounds, skin inflammation, or lumps.
- Persons with allergies or skin conditions such as dermatitis, etc.
- Do not use this device if you are pregnant or postpartum, as it may lead to skin concerns due to hormonal imbalances.
- Persons who suffer from seizures or epilepsy, congestive heart failure, or those undergoing radiation therapy or chemotherapy.

Do not use this device if your skin is sensitive to light. If you are taking photosensitizing agents or any other doctor-prescribed medication, please check the patient information leaflet before you start using this device. If the patient information leaflet states that it can cause photo-allergic or photo-toxic reactions, or if it specifies avoiding sun exposure when taking the medicine, please do not use the device. This device is not intended for use by persons (including children) with reduced physical, sensory, or mental

capabilities, or lack of experience and knowledge unless they have been given supervision, or instruction concerning the use of the device by a person responsible for their safety.

Cautious use of this device is advised if you have certain eye-related conditions, such as congenital retinal disorders.

If you have diabetes, please consult a doctor before using this device.

Discontinue use and consult a doctor if you experience any discomfort or abnormality while using the device. Prior to usage, always verify the condition of the device. Do not use it if it's damaged or broken.

⚠️ PRECAUTIONS

Use and store the device in a dry indoor environment.

Keep this device out of the reach of children.

Regularly clean the device to extend its lifespan.

Do not incinerate this device to prevent the risk of explosion.

When charging, be sure to use a standard USB adapter with the following specifications: DC 5V 1A. This helps prevent potential issues such as damage from high voltage or ineffective charging due to low current. This device has low battery voltage protection. If the charging indicator light flickers continuously and the device shuts down, the device battery is running low. Please make sure to fully charge the battery before using the device again.

If left unused for an extended period, it's recommended to charge the device once every six months.

Device Description TB-2386

Device Specifications

Name: GLAM LED Facial Mask

Model No.: TB-2386

Input: 5V 1A

Type of Battery: Li-polymer 3.7V 800mAh

Executive Standard: Q/XTM012-2020

Device Components (Picture 1)

① Main Unit
② Charging Cable
③ Instruction Manual

Picture 1

Device Features (Picture 2)

① Touch Switch
② Sensor Switch
③ Charging Indicator Light
④ LED Lamps (Light Therapy Mode):
- Red Mode, Yellow Mode, Blue Mode
⑤ USB Type-C Charging Port
⑥ Ear Hooks
⑦ Goggles

Picture 2

Preparing for Use TB-2386

Charging the Device

- Charging the device takes about 2.5 hours. A full battery allows for >60 minutes of continuous use.
- While charging, the charging indicator will flash a red light. Once fully charged, the charging indicator will light up continuously blue.
- This device has a low battery reminder. When the battery is low, the charging indicator light flashes blue 3 times with 3 'D' sounds. Please charge the device promptly. If the charging indicator light flashes red with 5 'D' sounds, the device will automatically shut down.
- This device does not come with an adapter. Please ensure that you use an adapter with an output of at least 1A. Turn ON/OFF
- Before starting, ensure that the device is sufficiently charged.
- Untie the ear hooks to activate the control switch and put the device into standby mode. Place the mask on your face to automatically start the light therapy mode.
- Remove the mask from your face to deactivate the light therapy mode and return to standby mode. If you have finished your session, fold the ear hooks to turn off the device. (Picture 3)

Picture 3

How to Use TB-2386

To ensure the correct operation of this device, please read the instructions carefully before use.

Before using the device for the first time, please remove the protective adhesive film from the front of the mask. (Picture 4)

Using the Device (Picture 5)

- Remove makeup and cleanse your skin.
- Unfold the ear hooks and place the mask on your face. Ensure the mask is properly in place by adjusting it to fit snugly around the eyes to prevent light leakage. The mask will automatically activate the light therapy mode.
- Press the touch switch twice to select the light therapy mode you want to use.
- Each mode operates in a 15-minute cycle. After 15 minutes, the device automatically deactivates the light mode and enters the standby mode.
- Remove the mask, fold the ear hooks, store the device, and proceed with your skincare routine.

※ Please close your eyes when placing or removing the device from your face to protect them from the LED lights.

Picture 4

Picture 5

Picture 6

Mode Switching: If the device is switched off without running out of battery, it will automatically activate the mode used in the last session when turned on again. However, if the device runs out of battery, it will revert to the default mode after being recharged. With the device activated in any mode, double-press the touch switch to cycle through red, blue, and yellow lights. (Picture 6)

Picture 6

Mode Function Description

※ You can choose the mode that best suits your skincare needs.

Mode	Firming Mode	Anti-Acne Mode	Brightening Mode
Technology	Red Light	Blue Light	Yellow Light
Effective Functions	Helps reduce wrinkles and fine lines. Helps stimulate the production of natural collagen to firm and plump the skin.	Helps reduce acne and prevent future breakouts by eliminating acne-causing bacteria. Gradually fades acne scars, providing a more even skin tone. Soothes skin redness and promotes a healthier complexion.	Soothing and anti-inflammatory, it helps reduce redness and inflammation. Helps heal and rejuvenate the skin, while also reducing sun damage. Helps lighten pigmentation. Prolonged usage can lead to a brighter and more even skin tone.
Use Frequency	For best results, use this device daily. We recommend selecting one mode per day and switching to another mode the following day. It is not advisable to repeat the same light therapy mode exposure.		

Product Management TB-2386

Safeguards

- Do not place or store this device where it can fall or in other situations that may cause external impact to the device.
- Store this device out of the reach of children.

Cleaning & Maintenance

- Ensure that the device is switched off before cleaning it.
- Do not rise or immerse this device in water or other liquid substances.
- Clean the surface of this device with a clean, slightly dampened cloth.
- Never use steel wool, abrasive cleaning products, chemicals, or cosmetic detergents to clean the device.

Storage

This device is equipped with a lithium battery. Keep it away from sunlight and avoid placing it near heat sources to prevent deformation, malfunction, and discoloration. Do not store it in humid or high-temperature environments. If you don't use the device frequently, ensure to charge it once every 6 months.

Disposal Information

Lithium batteries contain substances that can be harmful to the environment. To ensure responsible disposal and environmental protection: never incinerate this device, as it may pose the risk of explosion; do not dispose of this device with regular household waste. Instead, take it to an authorized recycling facility.

After-Sales Service TB-2386

Troubleshooting

Problem	Possible Cause	Solution
The device is not working properly.	The battery is running low or out of power.	Charge the device as guided in this instruction manual.
Other problems.		Please consult with your local distributor.

Note: If the device still does not work properly after using the above solutions, please stop using the device and contact the store or local distributor.

Warranty

※ Make sure that you receive an unopened box.

※ This warranty does not cover damages caused by the customer. In this case, the customer is responsible for covering the repair cost.

If you need to learn more about our products or if you need any information, please contact your local distributor or visit www.beauty.com.

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